## UV point irradiation device Laboratory-kit LKUVP-60 + UVPC

### Focal sizeΦ12 Easy UV Sterilization!



#### ♦ Feature ♦


- 1). Since it is a kit, user can easily use the UV irradiation device.
- 2). UV irradiation of 375 nm is possible.
- 3). UV irradiation of  $\Phi$ 13 can be easily reproduced.
- 4). The manual lifting platform allows for easy positioning.

# <u>Heat-tech</u>



[ Power:1500V-2W ]

It was adjusted to the position of 30 mm with a manual lifting platform. Rated point diameter can be visually confirmed.





Example of lab kit assembly.

(  $\$  The lab kit is delivered as individual components. )


[Lab kit packing product]


- 1 Power cable for heater controller
- (2) Power controller UVPC-1500V
- 3 Test stand
- (4) Heater mounting bracket
- (5) Ultraviolet point type irradiation device UVP-60/f30/1500v-2W/253.7nm With an output of 1500V-2W, it emits germicidal ultraviolet rays of 253.7nm.
- 6 Manual lift.


Up and down 80 mm by 16 rotating knobs, the accuracy is 0.5 mm. Manual lift. Platform surface is 100mm x 100mm,

O UV protection glasses

The irradiation point can be safely confirmed visually.







No-touch High Temperatures Hating

# **Heat-tech**

### Heat-tech Co., Ltd.

#### https://heater.heat-tech.biz/

International Medical Device Alliance IMDA 1–6–5 Minatojima Minamimachi Chuoku Kobe 650–0047 Japan TEL 81–78945–7894 FAX 81–78945–7895 E-mail info@heat-tech.biz